

Mesoscale light sheet microscopy for imaging cm-sized cleared tissue slices

Sharika Mohanan^{1*}, Steven Moreno¹, Callie Lorimer¹, Camilla Olianti², Eline Huethorst³, Erin Boland³, Leonardo Sacconi², Godfrey Smith³, Caroline Mullenbroich¹



¹ School of Physics and Astronomy, University of Glasgow, UK

² European Laboratory for Non-Linear Spectroscopy, Florence, Italy


³ School of Cardiovascular and Metabolic Health, University of Glasgow, UK

*sharika.mohanan@glasgow.ac.uk




1. Introduction


- The mesoSPIM (Mesoscale Selective Plane Illumination Microscope) is an open-source hardware light-sheet system [1].
- Built to image large cleared samples.
- Digitally scanned light-sheet with dual excitation arms and orthogonal detection arm.
- We characterise the mesoSPIM across its large field of view (FOV) and apply it to image cm-sized cleared rabbit tissue slices.




2. ASLM and zoom



3. Gaussian beam characterisation





4. Aberrations across FOV



5. Imaging of cleared tissue slices

- 500 μm tissue slices were excised from the left ventricle of the rabbit heart.
- Tissue clearing performed using hydrogel tissue transformation-based CLARITY protocol [4,5].
- CLARITY creates highly transparent tissue however with low structural integrity in the heart.
- We observe tissue expansion by × ~1.7.



- Cleared tissue mounted between two quartz slides immersed in refractive index-matching solution (EasyIndex).
- Sample angled at 45° with respect to detection axis.
- Tissue imaged in autofluorescence mode.
- Raw dataset deskewed to obtain XY projection.

6. Future work

- Improving the imaging quality of the system by using objectives that are aberration free across the entire 10.7 mm FOV at 2× magnification.
- Implement fluorescence staining on cleared tissue using Wheat Germ Agglutinin to demarcate myocardial cells.
- Assess and quantify structural differences in healthy vs. diseased heart tissue.

REFERENCES

- [1] Voigt *et al.*, Nature Methods (<https://doi.org/10.1038/s41592-019-0554-0>)
- [2] Deal *et al.*, Nature Protocols (<https://doi.org/10.1038/s41596-022-00706-6>)
- [3] Theer *et al.*, Nature Methods (<https://doi.org/10.1038/nmeth.3102>)
- [4] Tomer *et al.*, Nature Protocols (<https://doi.org/10.1038/nprot.2014.123>)
- [5] Olianti *et al.*, Prog Biophysics Mol Biology (<https://doi.org/10.1016/j.pbiomolbio.2021.07.012>)

FUNDING



- SM1

Cant change colormap values as
it is automatically set on psfj.

Sharika Mohanan,
2023-11-05T13:28:02.178
- SM2

Placeholder image. Will change
to correct image tomorrow.

Sharika Mohanan,
2023-11-05T13:29:07.845