# Imaging pipeline for cm-sized tissue slices using mesoscale light sheet microscopy

Dr Sharika Mohanan School of Physics and Astronomy University of Glasgow


(sharika.mohanan@glasgow.ac.uk)



## Motivation – Structural imaging of cardiac tissue

- Structural differences in healthy and diseased rabbit hearts
- Imaging tissue slices using protocols such as CLARITY
- Lateral extent of tissue much larger than thickness

Simplify imaging protocol






## Cardiac tissue – sample preparation

New Zealand
white rabbit
Left Ventricle

Vibratome
500 µm – 3 mm

Tissue
clearing
labelling



- Tissue expansion by 1.6x
- Tissue labelling
   Wheat Germ Agglutinin Alexa 488

Steven Moreno

Dr. Eline Huethorst

Dr. Erin Boland

Dr. Camilla Olianti












## Light sheet microscopy



#### Advantages:


- Widefield imaging
- Reduced photobleaching
- Large samples

#### **Considerations:**

- Sample mounting
- Data Handling



## Mesoscale Selective Plane Illumination Microscopy (MesoSPIM)



#### **Excitation arm:**

- Axially scanned light sheet (ASLM)
- Digitally scanned light sheet
- Excitation NA: 0.1


#### Detection arm:

- MVX10 Olympus zoom body
- Photometrics Kinetix camera (29.4 mm diagonal FOV)

Magnification	FOV (mm)
2x	14.7
4x	7.35
6.3 x	4.6




## MesoSPIM






## Tissue slice imaging – sample mounting





## Tissue slice imaging – sample mounting






Steven Moreno


- Mounted on a x,y,z,θ stage
- Tissue placed between two quartz slide "sandwich"
- Outer and inner "sandwich" cuvette filled with refractive index matching solution



## Tissue slice imaging





## **Oblique Scanning**






# Tissue slice imaging – oblique scanning








Scale bar = 500 µm

# Tissue slice imaging – oblique scanning









Scale bar = 500 μm FOV: 9000 x 3000 x 300 μm

## Tissue slice imaging – oblique scanning





#### Computational time for shearing

Lateral shift algorithm in MATLAB

Data size: 11 GB

3200x3200x500 pixels

Time:  $312 \pm 79 \text{ s}$ 


#### Spatial resolution in sheared images

	xy FWHM, μm	z FWHM, µm
3D z-scan	5.07 ± 0.76	6.51 ± 0.46 μm
Sheared slice scan	$7.22 \pm 0.43$	6.27 ± 1.08 μm

#### Mechanically Sheared Axially Swept Light-Sheet Microscopy

JINLONG LIN<sup>1,2</sup>, DUSHYANT MEHRA<sup>1,2</sup>, ZACH MARIN<sup>1,2</sup>, XIAODING WANG<sup>1,2</sup>, HAZEL M. BORGES<sup>1,2</sup>, QIONGHUA SHEN<sup>1,2</sup>, SEWERYN GAŁECKI<sup>1,2,3</sup>, JOHN HAUG<sup>1,2</sup>, AND KEVIN M. DEAN<sup>1,2,\*</sup>.


Lin et al, Bioarxiv, 2024




<sup>\*</sup>kevin.dean@utsouthwestern.edu




## Cardiac Structural Imaging



Scale bar =  $500 \mu m$ 



Scale bar = 250  $\mu$ m



Scale bar =  $25 \mu m$ 



### Conclusion

#### Advantages of oblique scanning:

- Reduced computational overhead
- Reduced tiling
- Maximum lateral extent (x') of tissue is determined by the working distance of objective
- No extra hardware or optics required for implementation



## Acknowledgements

#### University of Glasgow

#### Physics and Astronomy:

Dr. Caroline Müllenbroich Steven Moreno Giedre Astrauskaite

Dr. Ryo Kinegawa

Lewis Williamson

#### Institute of Clinical Physiology – CNR, Florence, Italy

Dr. Leonardo Sacconi

#### European Laboratory for Non-Linear Spectroscopy, Florence, Italy

Dr. Camilla Olianti

#### Cardiovascular and Metabolic Health:

Prof. Godfrey Smith

Dr. Eline Huethorst

Dr. Erin Boland














## MesoSPIM characterisation

#### Field of View






#### Resolution





Lateral resolution:  $5.07 \pm 0.76 \,\mu\text{m}$ Axial resolution:  $6.51 \pm 0.46 \,\mu\text{m}$ Mean of 70 beads across FOV





## Cardiac tissue – sample preparation





Riddell et al, Cardiovascular research, 2020





Dr. Eline








## Cardiac tissue – sample preparation

New Zealand
white rabbit
Left Ventricle

Vibratome
500 µm – 3 mm

Tissue
clearing

Iabelling







- Langendorff perfusion
- Tissue fixation in PFA
- Embedded in agar

Steven Moreno

Dr. Eline Huethorst

Dr. Erin Boland

Dr. Camilla Olianti











## Header

