

A pipeline of clearing, mounting, and mesoscale imaging of sliced cardiac tissue structure




Steven M Moreno*, Sharika Mohanan, Callie Lorimer¹, Camilla Olianti², Eline Huethorst³, Erin Boland³, Leonardo Sacconi², Godfrey Smith³, Caroline Müllenbroich¹

¹ School of Physics and Astronomy, University of Glasgow, UK

² European Laboratory for Non-Linear Spectroscopy, Florence, Italy

³ School of Cardiovascular and Metabolic Health, University of Glasgow, UK

*s.moreno.1@research.gla.ac.uk



1. Introduction to Optical System

- Observation of axially deep cellular tissues limited by the scattering of light [2].
- Optical clearing renders mesoscale tissues (mm³ to cm³) transparent, amenable to optical imaging.
- Light-sheet microscopy (LSM) permits isotropic resolution, large field of views, high throughputs of examined tissues [2].
- Preparation, clearing, staining, and mounting of the delicate cardiac tissue slices difficult to repeat.
- Reliable imaging results requires careful, reproducible protocols to mitigate tissue damage/alteration in vitro compared to in vivo.





Figure 1: mesoSPIM version 5 LSM at the University of Glasgow, Glasgow, United Kingdom.

2. Optical Pathways of mesoSPIM




mesoSPIM Excitation Arm (Figure 2a):

- Focus of the Gaussian beam shifted across sample using ETLs using Axially Swept Light-sheet Microscopy mode (ASLM, figure 2b) [2].
- ETL sweep synchronised with camera Photometrics Kinetix sCMOS camera (diagonal = 29.4 mm) rolling shutter feature.


mesoSPIM Detection Arm (Figure 2c):

- MVX-10 Olympus zoom body changes the FOV from 20.8 mm (1 \times) to 3.3 mm (6.3 \times) [2].
- MVXPLAP01 Olympus Objective focuses onto illuminated sample with 0.65 NA, 65mm WD




3. Cardiac Tissue Fixation and Slicing

Tissue Extraction and Fixation



- Extracted cardiac samples undergo Langendorff perfusion, expelling remaining blood in the organ (Figure 3a)
- Cardiac Fixation done using Paraformaldehyde (PFA), stored in Phosphate Buffered Saline (PBS) (Figure 3b) [1]

Tissue Slicing



- Fixated tissue is embedded in agar block
- Block glued onto base platform of Vibratome tissue slicer (Figure 3c)
- Tissue in agar sliced at desired thickness (400-3000 μm) (Figure 3d)
- Excess agar removed; slices returned to PBS storage.

4. Tissue Clearing Protocols: CLARITY

CLARITY Protocol [1,3]

(Clear Lipid-exchanged Acrylamide-hybridized Rigid Imaging Tissue hYdrogel)

- Washed in PBS to remove non-tissue molecules (PFA, blood, bacterium, etc.)
- Wash in Hydrogel solution (paraformaldehyde and hydrogel monomers), ensuring homogeneity.
- Degassing in N₂ polymerizes monomers into a hydrogel mesh.
- Wash in clearing solution (detergent with weak acid) to expel lipids.
- Changing out of clearing solution every couple of days (for up to 6 months)




Figure 4a: 400 μm Cardiac Tissue Sample Before and After CLARITY Tissue Clearing Protocol (± 20 weeks)





Figure 4b: CLARITY Tissue Clearing Protocol Steps and Timeline. [1,3]

5. Staining and Mounting of Cleared Samples



Tissue Staining

- Slice cardiac samples stained using Wheat Germ Agglutinin Conjugates (WGA-Fluorescein, Alexa Flour 488) (Figure 5a)
- Washed in 1x Phosphate Buffered Saline, Triton-X100 solution [1]
- Concealed from light throughout staining to prevent photobleaching (Figure 5b)

Tissue Mounting

- Cleared tissue mounted between two quartz slides immersed in Refractive Index matching solution (RI = 1.45)
- Sample angled at 45° from excitation/detection paths
- Angled position prevent frame from blocking optical pathways




Figure 5c: 1000 μm Cardiac Tissue Sample in custom mount before immersion in RI matching solution, imaging in mesoSPIM




Figure 5d: 3D model of mount at 0° [A] and 45° Orientation [B].

6. Imaging of Cleared Tissue Slices

- Z-stacks of stained samples recorded and processed
- Sliced with thickness of 500 μm , 1mm, 3mm produced high resolution images
- De-skewing program applied to obtain images stacks at 0°




Figure 6a: De-skewed Image of 500 μm Healthy Cardiac Tissue, x6.3 mag, WGA-Fluorescein Staining
Figure 6b: Mounted Slice Imaged in mesoSPIM

7. Future work

- Modify existing Detection Pathway using new latest open-source design to mitigate aberrations, vignetting associated with expanded FOV
- Implement staining on samples from diseased hearts to assess structural changes relative to healthy hearts
- Proceed with quantitative, structural analysis of cardiac tissue for images obtained.
- Examine clearing results from alternative clearing protocols (SHIELD, CUBIC-L/RA), assessing image quality relative to CLARITY

REFERENCES

- Olianti, Camilla, et al. *Progress in Biophysics and Molecular Biology* 168 (2022): 10-17.
- Voigt, Fabian F., et al. *Nature methods* 16.11 (2019): 1105-1108.
- Tomer et al, *Nature protocols* 9.7, (2014): 1682-1697.

FUNDING

